



Pearl II™



SOS Button

Safe Zones

Fall Detection

Location Updates

Two-way Voice Calls

Audio Reminders

Easy Charging

Auto Answer

Introducing Oysta Pearl II™

The Pearl II has been designed from direct feedback from Vulnerable Independent People (VIPs) and their care circles.

It retains many of the features of the proven and loved Pearl+ such as an SOS button, fall detection, two-way voice calls and audio reminders. It has a new streamlined design and charging cradle to avoid fiddly wires for those with hand mobility issues.

The Pearl II is a mobile digital solution which enables early intervention, preventative and self-care for VIPs who want to remain independent.

The Oysta Pearl II connects to our care platform, IntelliCare™, so the care circle will be able to set safety parameters, locate VIPs and if needed, guide them home or direct help to their location.

The Pearl II is ideal for:

- Active VIPs that have started to experience memory loss/disorientation and may need locating
- VIPs who would benefit from a simplistic device design
- VIPs with visual impairments; the Pearl II has only three main buttons and with braille on them
- VIPs that lack dexterity to use standard chargers

Benefits for Vulnerable Independent People (VIPs)

One Touch SOS Button: An easy to use button when the VIP is in distress or simply needs peace of mind. The Pearl II connects the VIP to their chosen emergency contact or Alarm Receiving Centre (ARC) and makes a hands-free, two-way voice call.

In the event of a fall an SOS call is automatically activated.

Two-way Voice Calls: The Pearl II has a loudspeaker and sensitive microphone to allow clear voice communication.


Audio Reminders when the device has a low battery to prompt the VIP to charge the device.

Cradle Charging means that the VIP only needs to drop it into the charging cradle to charge the device.


The **Auto Answer** option means that the call is automatically answered so there is no need to press a button which saves time in an emergency or a welfare check situation.




Alert raised
on Pearl II



Alert answered by
Alarm Receiving
Centre (ARC)



ARC receives
information
and location



Help is
sent to
the VIP

Benefits for the Care Circle

Fall Detection

Should the VIP have a fall it will automatically connect them to their chosen contact or Alarm Receiving Centre (ARC).

Location Updates

You can see on IntelliCare™ where your VIP is and if they need help, you can locate them and guide help to their location.

Safe Zones

Safe Zones are areas where the VIP is likely to travel. Should the VIP travel outside those safe zones anytime, or during specific times, of the day or night, a care alert is triggered.

WeatherWatch

A significant step towards preventative care; WeatherWatch notifies the care circle if the VIP is outside after dark or in extreme weather conditions. You can conduct a welfare check before problems occur, and if needed, guide them home using the location feature on IntelliCare.

Reminder Messages

Messages can be pre-set, flashing up on the Pearl II display screen to remind the VIP to take medication.

Whitelist of Contacts

Only authorised callers from the care circle can contact the Pearl II which eliminates spam calls to the VIP. The whitelist can be changed at any time to add/remove contacts.

Status Alerts

Through IntelliCare, you'll be notified if your VIP's device is switched on/off, has a low battery or if there has been 'no activity', indicating that the device has not been worn or used.

Intelligent Care Technology

Oysta helps Vulnerable Independent People (VIPs) stay safe at home and away from home. Our technology encourages VIPs to pursue an active lifestyle with a safety net in place.

The Oysta Pearl II™ is one of our telecare devices which protect VIPs when away from their home. We have also developed a range of SafeHome™ sensors to keep VIPs safe when they are at home and enable preventative care efforts of the care circle.

Our devices and sensors link into our easy-to-use care platform, IntelliCare™. IntelliCare is used by the care circle so that they have regular, automatic updates on activity and location enable trends and analysis to support effective care decisions.

Oysta's IntelliCare platform simplifies technology enabled care and transforms reactive care to personalised, preventative care, helping customers to deliver better outcomes.

SPECS	
Size	84mm x W 45mm x D 17mm
Weight	65g
FEATURES	
IP Rating	IP44
Speaker	Yes
Microphone	Yes
TECH	
Signal	2G/3G/4G
GSM Freq.	Quad-band 850/900/1800/1900 MHz
GPS	GPS/AGPS/LBS/Wifi
Battery	Li-Polymer 850mAh
Operating temp	-20°C ~ +55°C
Display	250mm x 278mm
Sensors	Accelerometer

Call us on
01295 530 101
to book a live
demonstration of
the Pearl II™ and
IntelliCare™



OYSTA

OYSTA TECHNOLOGY

Home Ground Barn, Pury Hill Business Park, Towcester, NN12 7 LS, UK.

Tel +44 (0) 1295 530 101

Email info@oysta-technology.com

www.oysta-technology.com



Version: Jul2022